



Virusolve[®]+ EDS

High Level Disinfectant / Cold Sterilant for Endoscopes & Medical Devices

Virusolve+ EDS is a superior, patented product for the decontamination of medical, surgical, dental and endoscopic instruments.

Primary features attributable to Virusolve+ EDS are:

- High Level Disinfectant / Cold Sterilant.
- No instrument damage.
- Excellent material compatibility.
- Safe.
- Odourless.
- Non-fixative.
- Non-corrosive.
- 5 minute soak at 20 °C.
- No activation required.
- Cost effective.
- Tested to EN14885 for disinfectant properties



Virusolve[®]+ EDS

A High Level Disinfectant / Sterilant Fluid for Endoscopes & Medical Instruments

General

Virusolve+ EDS is a superior high level disinfectant (HLD) (Sterilant) that is designed for the decontamination of endoscopic, surgical and dental instruments & medical devices.

Virusolve+ EDS is a patented proven HLD / Sterilant. The formulation does not contain any hazardous aldehydes or chlorine generating components. It is highly effective as a HLD / Sterilant against bacterial spores, mycobacteria, bacteria, fungi and viruses. The product is designed to eliminate the risk of crossinfection from spores, mycobacteria, viruses, fungi and bacteria.

Virusolve+ EDS is safe to use on ALL flexible and rigid endoscopes.

Virusolve+ EDS Microbiological Properties

| Active Against | Standards | Contact Time |
|----------------|-------------------------------------------------------------------------------------------------------------------------|------------------------------------------------|
| Spores | EN13704 EN14347 | 2 Minutes 1 Minute |
| Mycobacteria | EN14348 Mycobacterium terrae Mycobacterium tuberculosis EN14563 Mycobacterium terrae Mycobacterium avium | 1 Minute 1 Minute 5 Minutes 5 Minutes |
| Bacteria | EN1276, EN13727, EN14561 | 1 Minute 5 Minutes |
| Viruses | EN 14476 Polio & Adeno HIV, Hep B & C | 5 Minutes 1 Minute |
| Fungi / Yeasts | EN1275, EN1650 EN14562, EN13624 | 1 Minute 5 Minutes |

Packaging

Concentrate

- 125ml container
- 250ml container
- 1 L container
- 5 L container
- 10 L container
- 25 L container

Ready to Use

- 5 L container



Application

REPROCESSING OF INSTRUMENTS - HIGH LEVEL DISINFECTION


For ALL critical and semi-critical instruments / medical devices including heat sensitive instruments:

For Manual Disinfection of instruments and equipment:


Stage 1:

- Make sure that instrument has been thoroughly cleaned (with Viruzyme family or Virudet N Manual) and rinsed.

Stage 2:




Make Virusolve+ EDS solution at 25ml/L (2.5%) with water.



Soak for 5 minutes making sure that all internal channels are flushed through. This ensures sporicidal, mycobactericidal, virucidal, bactericidal, and fungicidal performance is attained.

Note: Ensure that all spaces and recesses are completely filled with Virusolve+ EDS fluid and no air pockets exist



Rinse in sterile water (repeat) (2 rinses are recommended) and dry with sterile material.

The medical device is now ready for use

Re-usage for Disinfection

Virusolve+ EDS solution has demonstrated efficacy in the presence of organic soil contamination and microbiological burden during reuse. The active constituents' concentration of Virusolve+ EDS solution during its use-life must be verified by the Virusolve+ EDS Concentration Test Kit or Test Strips.

Virusolve+ EDS must be discarded after 14 days even if the active constituents are at or above the MEC. For MEC please refer to PIBS 08-006.



Distributed by:

Amity International,
Libra House, West Street,
Worsborough Dale,
Barnsley, S70 5PG, UK.

Tel: +44 (0) 1226 770787
Email: sales@amityinternational.com
www.amityinternational.com

© Amity Limited all rights reserved (dec 2014)

